FUJIFILM**Wako**Code No. 293-77601 (10 purifications)
299-77603 (2 purifications)

For Genetic Research

MagCapture™ Exosome Isolation Kit PS**【Introduction】**

Exosomes and other extracellular vesicles (EVs) are small membrane vesicles containing protein, mRNA, microRNA, DNA, and lipids, which are secreted by various cells and are stable in body fluids including blood, saliva, urine, cerebrospinal fluid (CSF) and breast milk. These EVs have been recognized as messengers of cell-to-cell communication, and biomarkers for various diseases. Conventionally, ultracentrifugation, affinity purification using antibody to surface antigen, and precipitation with polymer reagents are used for isolation of exosomes. However, these methods are not fully satisfying in all performances criteria such as recovery efficiency, purity and operability.

MagCapture™ Exosome Isolation Kit PS adopts a novel affinity purification method using magnetic beads and phosphatidylserine (PS)-binding protein (PS affinity method). This kit can easily isolate high purity exosomes and other EVs from cell culture medium and body fluids at high yield by a normal microcentrifuge. If higher purity exosomes are needed, please use the 10,000 × g supernatant as sample. This kit enables the isolation of exosomes and other EVs as intact forms because the captured EVs are eluted from magnetic beads with metal-chelating reagent at neutral pH. The isolated intact exosomes and other EVs can be used for various applications including electron microscopic analysis, nanoparticle tracking analysis, administration of EVs and analysis of molecular constituents such as proteins, lipids, or nucleic acids.

【Features】

- High purity exosomes and other EVs are isolated by a novel affinity purification method (PS affinity method).
- The purity and yield of exosomes are higher than ultracentrifugation.
- This kit enables the purification of exosomes and other EVs as intact forms, and can be used for various applications.
- Exosomes can be isolated with magnetic beads.
- No ultracentrifugation is required.

【Kit contents】

This kit includes 7 components.

· 2 purifications

- | | |
|---|------------------|
| (1) Streptavidin Magnetic Beads | 120 μL × 1 tube |
| (2) Biotin-labeled Exosome Capture | 20 μL × 1 tube |
| (3) Exosome Capture Immobilizing Buffer | 7 mL × 1 bottle |
| (4) Exosome Binding Enhancer (× 500) | 100 μL × 1 tube |
| (5) Washing Buffer | 30 mL × 1 bottle |
| (6) Exosome Elution Buffer | 1 mL × 1 bottle |
| (7) Reaction Tubes | 4 tubes |

- 1/24 -

· 10 purifications

- | | |
|---|-------------------|
| (1) Streptavidin Magnetic Beads | 600 μL × 1 tube |
| (2) Biotin-labeled Exosome Capture | 100 μL × 1 tube |
| (3) Exosome Capture Immobilizing Buffer | 35 mL × 1 bottle |
| (4) Exosome Binding Enhancer (× 500) | 500 μL × 1 tube |
| (5) Washing Buffer | 75 mL × 2 bottles |
| (6) Exosome Elution Buffer | 5 mL × 1 bottle |
| (7) Reaction Tubes | 22 tubes |

【Storage】

Store at 2-10 °C

【Additional required materials】**Reagents :**

- 1) TBS (As necessary)

Equipment :

- 1) Microcentrifuge (Max > 10,000 × g)
- 2) Vortex mixer
- 3) Tabletop centrifuge
- 4) Magnetic stand
- 5) Rotator or tube mixer
- 6) Micro pipette
- 7) Pipette tip
- 8) Microcentrifuge tube (1.5 mL)
- 9) Centrifuge tube (15 mL) (As necessary)
- 10) Centrifuge tube (50 mL) (As necessary)
- 11) Ultrafiltration unit (Vivaspin6, M. W. cut off : 100K, Code No. VS0641 or Vivaspin20, M. W. cut off : 100K, Code No. VS2041) (As necessary)

【Precaution for use】**1. Equipment**

Please use sterile or DNase and RNase free microcentrifuge tubes and pipette tips. We recommend the use of gloves and a mask to avoid contamination of DNase and RNase.

2. Reagents


If you analyze the nucleic acid (RNA, DNA) after isolation of exosomes and other EVs, please handle reagents carefully and keep them as sterile as possible. Additionally, if you analyze the mass spectrometry of proteins after isolation of exosomes and other EVs, please be careful of contamination from other proteins.

3. Handling of biohazardous waste

Please treat human serum, human plasma, and human tissue samples as infectious samples. Biohazardous wastes, which include waste liquid and equipment such as microcentrifuge tubes, pipette tips, and gloves, must be disposed of according to the guidelines of the institution.

- 2/24 -

【Outline of procedure】

【Procedure】

1. Preparation of sample

This is the section to prepare samples. When exosomes and other large EVs are needed, prepare **1,200 × g supernatant** as a sample. Additionally, when highly purified exosomes are needed, prepare **10,000 × g supernatant** as a sample^{※1}. This protocol for sample preparation is set for cell culture medium, serum, and plasma. When other body fluids are used, please examine the appropriate pretreatment protocol by referring to the protocol for serum and plasma.

※1 To remove Large EVs such as apoptotic vesicles and microvesicles from a sample, use of a centrifugal filter unit (Millipore, Ultrafree-MC, GV 0.22 μm sterile, Product Code : UFC30GV0S) on the “supernatant” obtained from the sample by centrifugation at 10,000 × g for 30 minutes is recommended. The filtrate can be used as a test sample.

In the case of cell culture medium

- 1) Culture the cells under the appropriate condition^{※2}.
- 2) Harvest cell culture medium.
- 3) Centrifuge the cell culture medium for 5 minutes at 300 × g and 4 °C to remove cells and debris.
- 4) Transfer the supernatant from step 3) into a new tube.
- 5) Centrifuge 4) at 1,200 × g for 20 minutes at 4 °C to remove the cell debris.
- 6) Transfer the supernatant from step 5) into a new tube. (**1,200 × g supernatant**)
- 7) Centrifuge 6) at 10,000 × g for 30 minutes at 4 °C to remove the large EVs.
- 8) Transfer the supernatant from step 7) into a new tube. (**10,000 × g supernatant**)

【Optional: Concentration of cell culture supernatant by ultrafiltration and cautions】

If the volume of **1,200 × g supernatant** or **10,000 × g supernatant** is over 1 mL, concentrating the sample by using ultrafiltration unit (Vivaspin6, M. W. cut off : 100K, Code No. VS0641 or Vivaspin20, M. W. cut off : 100K, Code No. VS2041)

is recommended. Please concentrate until the sample is less than or equal to 1 mL in accordance with the manufacturer’s instruction manual. The maximum sample volume of ultrafiltration treatment is 50 mL (50-fold concentration). Exosomes and other EVs can be efficiently recovered from concentrated sample.

Cautions for ultrafiltration concentration

Since extracellular vesicles may be adsorbed to the centrifugal ultrafiltration unit and the amount recovered may be reduced, add EV-Save™ Extracellular Vesicle Blocking Reagent (Code : 058-09261) (sold separately) to the sample for preventing loss of them. However, EV-Save™ contains polymers, its use is not recommended for proteomic analysis samples.

※2 Please examine appropriate culture conditions depending on cell lines. Since the amount of extracellular vesicles in the medium is small, please increase the culture scale as much as possible. (Example: 20 mL or more is better) In addition, to ensure that isolated exosomes which originate from your cells, please use exosome-depleted FBS. Exosome-depleted FBS is prepared by ultracentrifugation (Example: at 110,000 × g for 18 hours) or is commercially available.

In the case of Serum and Heparin Plasma

Please centrifuge serum and heparin plasma used in this section at 1,200 × g beforehand or prepare the processed sample.

When exosomes and Large EVs are needed, 1,200 × g supernatant is used as a sample for purification. In that case, a floating matter may be observed in 1,200 × g centrifuged samples. When it was observed, centrifuge again at 1,200 × g and use it as a sample for exosome and Large EVs purification.

- 1) Centrifuge a sample at 10,000 × g for 30 minutes at 4 °C to remove the large EVs^{※3}.
- 2) Transfer the supernatant from step 1) into a new tube. … (**10,000 × g supernatant**)

※3 When Large EVs are needed, use ppt of Large EVs obtained by centrifugation at 10,000 × g as sample after suspending it with 500 μL to 1 mL of TBS.

In the case of EDTA Plasma and Citrated Plasma


Please centrifuge EDTA plasma and citrated plasma used in this section at 1,200 × g beforehand or prepare the processed sample.

In the case of EDTA and citrated plasma, the anticoagulant contained in the sample inhibits the binding of extracellular vesicles to Exosome Capture Beads. Therefore, adjust the heparin sodium solution just before the affinity reaction, add it to the centrifuged sample, and proceed to the purification step. Details are described in the section on **EDTA Plasma and Citrated Plasma** in 3. Affinity Reaction.


- 1) Centrifuge a sample at 10,000 × g for 30 minutes at 4 °C to remove the large EVs^{※3}.
- 2) Transfer the supernatant from step 1) into a new tube. … (**10,000 × g supernatant**)

[Flowchart of Sample Preparation (Section 1)]


Cell Culture Medium


Concentration (optional)^{※4}

Serum and Heparin Plasma

EDTA Plasma and Citrated Plasma

- ※1 When **exosomes and large EVs** are needed, use 1,200 x g sup. fraction as sample.
- ※2 When **exosomes** are needed, use 10,000 x g sup. fraction as sample.
- ※3 When **Large EVs** are needed, use ppt of Large EVs obtained by centrifugation at 10,000 x g as sample after suspending it with TBS.
- ※4 The concentration step is an option when using a large volume (~ 50 mL) of cell culture supernatant as a sample for purification. However, since recovery efficiency improves, please perform it as much as possible.

2. Immobilization of Exosome Capture


This is the section to immobilize Biotin-labeled Exosome Capture reagent to Streptavidin Magnetic Beads. Please use the 1.5 mL Reaction Tube (included) in the downstream experiment which requires magnetic beads operation^{※4}.

- 1) Transfer 60 μL ^{※5} of Streptavidin Magnetic Beads into a new 1.5 mL Reaction Tube.
- 2) Add 500 μL of Exosome Capture Immobilizing Buffer into the 1.5 mL Reaction Tube from step 1) and suspend it by vortexing.
- 3) Spin down 2), then places it on the magnetic stand for 1 minute. Remove the supernatant after separation of magnetic beads and supernatant.

- 4) Add 500 μL of Exosome Capture Immobilizing Buffer and 10 μL of Biotin-labeled Exosome Capture into the 1.5 mL Reaction Tube from step 3) and mix by vortexing.
- 5) Mix for 10 minutes at 2-10 $^{\circ}\text{C}$ with rotator or tube mixer.
- 6) Spin down 5), then places it on the magnetic stand for 1 minute. Remove the supernatant after separation of magnetic beads and supernatant.
- 7) Add 500 μL of Exosome Capture Immobilizing Buffer into the 1.5 mL Reaction Tube from step 6) and mix by vortexing.
- 8) Spin down 7), then places it on the magnetic stand for 1 minute. Remove the supernatant after separation of magnetic beads and supernatant.
- 9) Repeat steps 7)-8) twice. **Exosome Capture-immobilized beads**

- ※4 Please only use the 1.5 mL Reaction Tube included in this kit because only these tubes can generate the high recovery rate of extracellular vesicles by preventing non-specific adsorption of magnetic beads to the tube wall.
- ※5 Please adequately suspend the magnetic beads by vortexing. Make sure homogeneous suspension of the magnetic beads is achieved before use.

[Flowchart of immobilization of Exosome Capture (Section 2)]

3. Affinity Reaction

This is the affinity reaction section of the Exosome Capture-immobilized beads with exosomes and other EVs in sample. This protocol is based on the use of 1.5 mL Reaction Tube (max volume : 1 mL). When exosomes and large EVs are needed, use **1,200 x g supernatant** as sample. Additionally, when highly purified exosomes are needed, use **10,000 x g supernatant** as sample. Please note that the affinity reaction protocol for EDTA plasma and citrated plasma are different.

In the case of cell culture medium, serum, and heparin plasma

- 1) Transfer sample (max volume : 1 mL)^{*6} to a new 1.5 mL microcentrifuge tube (not included in this kit) and add 1:500 volume of Exosome Binding Enhancer (× 500) to sample. Mix by vortexing. . . . **A sample including Exosome Binding Enhancer (× 500)**
- 2) Spin down **A sample including Exosome Binding Enhancer (× 500)**, then transfer it into the 1.5 mL Reaction Tube containing **Exosome Capture-immobilized beads**. Mix by vortexing.
- 3) Mix for 3 hours ~ at 2-10 °C with rotator.
- 4) Spin down 3), then places it on the magnetic stand for 1 minute. Remove the supernatant^{*7} after separation of magnetic beads and supernatant. . . . **EVs-binding beads**

In the case of EDTA/citrated plasma

- 1) Dissolve 10,000U of sodium heparin (Code : 085-00134) with purified water to make 1,000U/mL heparin sodium solution^{*8}.
 - 2) Transfer sample (max volume : 1 mL) ^{*6} to a new 1.5 mL microcentrifuge tube (not included in this kit) and add 1 : 200 volume of heparin sodium solution to sample. (final conc. : 5 U/mL)
 - 3) Add 1:100 volume of Exosome Binding Enhancer (× 500) to 2) and mix by vortexing. . . . **A sample including Exosome Binding Enhancer (× 500)**
 - 4) Spin down **A sample including Exosome Binding Enhancer (× 500)**, then transfer it into the 1.5 mL Reaction Tube containing **Exosome Capture-immobilized beads**. Mix by vortexing.
 - 5) Mix for 3 hours ~ at 2-10 °C with rotator.
 - 6) Spin down 5), then places it on the magnetic stand for 1 minute. Remove the supernatant^{*7} after separation of magnetic beads and supernatant. . . . **EVs-binding beads**
- ※6 Add the appropriate volume of TBS buffer to each sample to reach a volume of 0.5 mL to obtain the better mixture of the Exosome Capture-immobilized beads with sample. (Example : 100-200 μL → 500 μL)
- ※7 When further recovery of remaining extracellular vesicles is needed, please store the used supernatant.
- ※8 The maximum usage per kit (for 10 purification) is 100 μL. After preparation, please keep refrigerated.

4. Washing of EVs-binding beads

This is the washing section of **EVs-binding beads**. Please use Washing Buffer after adding Exosome Binding Enhancer (× 500) certainly^{*9}.

- 1) Add 1 : 500 volume of Exosome Binding Enhancer (× 500) to the required Washing Buffer^{*10} and mix by vortexing^{*11}.
- 2) Add 1 mL of Washing Buffer (+ Exosome Binding Enhancer (× 500)) into the 1.5 mL Reaction Tube containing **EVs-binding beads** and mix by vortexing.
In the case of EDTA/citrated plasma
After adding 1:200 volume of heparin sodium solution to 1 mL of Washing buffer (+ Exosome Binding Enhancer (× 500)) only at the first washing (final conc. : 5U/mL), suspend with a vortex mixer.
- 3) Spin down 2), then places it on the magnetic stand for 1 minute. Remove the supernatant after separation of magnetic beads and supernatant.

- 4) Add 1 mL of Washing Buffer (+ Exosome Binding Enhancer (× 500)) and mix by vortexing.
- 5) Spin down 4), then places it on the magnetic stand for 1 minute. Remove the supernatant after separation of magnetic beads and supernatant.
- 6) Repeat steps 4)-5) one more time.
- 7) Spin down again the 1.5 mL Reaction Tube from step 6), then places it on the magnetic stand for 1 minute. Remove the supernatant completely after separation of magnetic beads and supernatant^{*12}. . . . **Washed EVs-binding beads**

※9 If the addition of Exosome Binding Enhancer (× 500) to Washing buffer is omitted, captured exosomes and EVs will be eluted from the magnetic beads by Washing Buffer during washing section.

※10 The required volume of Washing Buffer is 3 mL per reaction.

※11 Please prepare Washing Buffer (+ Exosome Binding Enhancer (× 500)) just before use, because constituents may be precipitated after the addition of Exosome Binding Enhancer (× 500) overtime.

※12 Please remove Washing Buffer (+ Exosome Binding Enhancer (× 500)) in the tube completely, because remaining buffer causes reduction of elution efficiency of captured exosomes and other EVs.

5. Elution of extracellular vesicles

This step is to elution section of extracellular vesicles from **Washed EVs-binding beads**. By adding EV-Save™ Extracellular Vesicle Blocking Reagent (Code : 058-09261) (sold separately) to the Exosome Elution Buffer used in this section, the adsorption of extracellular vesicles to the tube can be suppressed and the recovery efficiency can be increased. However, EV-Save™ contains polymers, its use is not recommended for proteomic analysis samples.

- 1) Add 50 μL of Exosome Elution Buffer into the 1.5 mL Reaction Tube containing **Washed EVs-binding beads** and mix by vortexing^{*13}, then spin down.
- 2) Places it at room temperature for 10 minutes. After the reaction, resuspend by vortexing and spin down^{*13}.
- 3) Places the Reaction Tube from step 2) on the magnetic stand for 1 minute. Transfer the supernatant into a new 1.5 mL microcentrifuge tube (not included in this kit) after separation of magnetic beads and supernatant.
- 4) Add 50 μL of Exosome Elution Buffer into the 1.5 mL Reaction Tube from step 3) again, and mix by vortexing^{*13}, then spin down.
- 5) Places it at room temperature for 10 minutes. After the reaction, resuspend by vortexing and spin down^{*13}.
- 6) Places the Reaction Tube from step 5) on the magnetic stand for 1 minute. Transfer the supernatant into the 1.5 mL microcentrifuge tube from step 3) after separation of magnetic beads and supernatant, and then pool together^{*14, 15}. (Total: 100 μL)


※13 If the magnetic beads aggregate and can not be uniformly suspended with only a vortex mixer, suspend the beads evenly while applying tapping or pipetting operation. Please avoid excessive suspension using a tube mixer.

※14 **Exosome Capture-immobilized beads** are reusable after eluting EVs (up to 4 times). Therefore, when you need to recover more remaining EVs from samples, please repeat the section 3-5 to increase the recovered amount.

After recovering the eluate, please add the stored supernatant of section 3 (refer to ※7) into the 1.5 mL Reaction Tube containing **Exosome Capture-immobilized beads** from step 6) of section 5 and repeat the section 3-3) to the section 5 as necessary. When using EDTA/Citrated plasma, repeat from 3-5). This kit includes the required reagents to reuse a maximum of 4 times.

※15 There is a possibility that the pooled eluate have a small amount of magnetic beads. When extracellular vesicles are analyzed by Nanoparticle Tracking Analysis (NTA) or electron microscopic analysis, please filter pooled eluate by using 0.45 μm pore size filter unit (Example: Ultrafree[®] MC 0.45 μm, Code No. UFC30HV25, Millipore[®]) before analysis as needed.

[Flowchart of Affinity Purification (Section 3-5)]

[Optional]

1. RNA extraction from Purified Extracellular Vesicles

In the case of RNA purification from eluted extracellular vesicles, please add 100 μL of Exosome Elution Buffer to pooled eluate to prepare 200 μL of sample in total. Then, extract RNAs from extracellular vesicles with Wako's microRNA Extractor SP Kit (Code No. 295-71701) in accordance with the instruction manual.

2. Recycling of Exosome Capture-immobilized beads after eluting of extracellular vesicles

Used Exosome Capture-immobilized beads can be recycled for "Repeated extraction of extracellular vesicles from the

same sample (refer to ※14)" and "Purification of extracellular vesicles from culture supernatant sample of the same lot and body fluid sample of the same lot". Recycling is up to 4 times.

When purifying extracellular vesicles from samples of the same lot, prepare samples first and proceed from affinity reaction [section 3-1)].

When storing used Exosome Capture-immobilized beads, please add the Washing buffer in the kit or 1 × TBS containing 0.05 w/v% sodium azide to the Reaction Tube containing the used Exosome Capture immobilized beads and suspend it. After suspension, keep it at 2-10 °C.

As a record of our company, it has been confirmed that after 3 months elapsed Exosome Capture immobilized-beads can be reused. However, please use it as soon as possible.

[Related Products]

Code No.	Description	Size
297-79201	PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)	96 reactions
298-80601	PS Capture™ Exosome ELISA Kit (Streptavidin HRP)	96 reactions
297-79701	PS Capture™ Exosome Flow Cytometry Kit	300 reactions
290-80301	PS Capture™ Exosome Isolation Resin Kit	1 Kit (0.5 mL Slurry)
016-27061	Anti CD63, Monoclonal Antibody (3-13)	20 μL
012-27063		100 μL
018-27641	Anti CD63, Monoclonal Antibody (3-13),	25 tests
014-27643	Fluorescein Conjugated	100 tests
013-27711	Anti CD63, Monoclonal Antibody (3-13),	20 μL
019-27713	Biotin Conjugated	100 μL
011-27751	Anti CD63, Monoclonal Antibody (3-13),	25 tests
017-27753	Red Fluorochrome (635) Conjugated	100 tests
018-27761	Anti CD9, Monoclonal Antibody (1K)	20 μL
014-27763		100 μL
015-27771	Anti CD81, Monoclonal Antibody (17B1)	20 μL
011-27773		100 μL
058-09261	EV-Save™ Extracellular Vesicle Blocking Reagent	1 mL
290-35591	Magnet Stand	1 each
295-71701	microRNA Extractor SP Kit	50 reactions
317-90175	10 × TBS (pH 7.4)	500 mL

[References]

- 1) W. Nakai, *et al.*, *Sci. Rep.*, **6**, 33935 (2016).
- 2) A. Kraft, W81 Award. (2017).
- 3) S. Osawa, *et al.*, *Biochem. Biophys. Res. Commun.*, **488**, 232-238 (2017).
- 4) S. Nagashima, *et al.*, *J. Virol.*, **91**(22), e00822-17 (2017).
- 5) T. Yoshida, *et al.*, *Curr. Protoc. Cell Biol.*, **77**(1), 3.45.1-3.45.18 (2017).
- 6) S. Saito, *et al.*, *Sci. Rep.*, **8**, 3997 (2018).

- 7) H. Kobayashi, *et al.*, *Nagoya J. Med. Sci.*, **80**(2), 141-153 (2018).
- 8) Y. Obata, *et al.*, *JCI Insight*, **3**(8), e99680 (2018).
- 9) R.C. Middleton, *et al.*, *J. Extracell. Vesicles*, **7**(1), 1456888 (2018).
- 10) H. Kawahara, and R. Hanayama, *Biol. Pharm. Bull.*, **41**, 1119-1125 (2018).
- 11) M. Santiana, *et al.*, *Cell Host & Microbe*, **24**, 208-220 (2018).
- 12) T. Antes, *et al.*, *J. Nanobiotechnol.*, **16**, 61 (2018).
- 13) H. Kimura, *et al.*, *Clinical Cancer Research*, **25**(6), 1936-1947 (2019).
- 14) G. Gorji-Bahri, *et al.*, *Current Gene Therapy*, **18**, 336-350 (2018).
- 15) H. Xu, *et al.*, *Anal. Chem.*, **90**(22), 13451-13458 (2018).
- 16) H. Kuroda, *et al.*, *Mol Pharm.*, **16**, 292-304 (2019).
- 17) J. Juan, *et al.*, *Stem. Cell Res. Ther.*, **10**(1), 95 (2019).
- 18) L. Zhi, *et al.*, *Cell Death Discovery*, **5**, 79 (2019).
- 19) G.K. Patel *et al.*, *Sci Rep.*, **9**(1), 5335 (2019).
- 20) Y.L. Tai *et al.*, *J. Biomed. Sci.* **26**(1), 35 (2019).
- 21) Y. Sugihara *et al.*, *PLoS One.*, **14**(5), e0217394 (2019).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : +81-6-6203-3741
Facsimile : +81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : +49-2131-3111-0
Facsimile : +49-2131-311100
<http://www.wako-chemicals.de>

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MagCapture™ Exosome Isolation Kit PS

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